

Preparation of Chromatin Fibers from Tissue Culture Cells

This protocol requires a cytocentrifuge that dispenses unfixed cell or nuclei suspensions onto glass slides. (I use the Cytospin 4 from Shandon, Inc.)

1. Transfer cells from flask into a 50 ml tube and bring final vol to 40 mLs with 1X PBS.
2. Count cells using a hemacytometer.
3. Resuspend in hypotonic (75mM KCl or 0.8% Na Citrate) to $0.8-2 \times 10^5$ cells/ml and centrifuge onto charged (Superfrost plus) slides. Use double funnels; 4 min at 800 rpm, high acceleration. For *Drosophila* cells, do not use superfrost plus (charged) slides; use regular glass slides.
4.
 - 2×10^5 cells/ml for *drosophila* S2 or Kc cells - single chamber
 - 3.5×10^5 cells/ml for *drosophila* S2 or Kc cells - double chamber
 - $0.5-1 \times 10^5$ cells/ml for human lymphoblasts - single chambers
 - $1-1.5 \times 10^5$ cells/ml for human lymphoblasts - double chambers
 - 0.5×10^4 cells/ml for mammalian primary fibroblasts - single chambers
 - 1×10^4 cells/ml for mammalian primary fibroblasts - double chambers
 - 1×10^4 cells/ml for mammalian immortalized fibroblasts - single chambers
 - 2×10^4 cells/ml for mammalian immortalized fibroblasts - double chambers
5. After cytopinning, immerse slides in a coplin jar of Lysis Buffer for 10-15 minutes at RT. You may have to vary the incubation time depending on the cell line (start with 8 minutes; if fibers are not stretched enough, increase time in Lysis Buffer. Alternatively, you can keep the shorter time (8m) but increase NaCl concentration to 625mM – be aware that you may lose proteins from the chromatin with increased salt.)
6. Slowly lift slide from lysis buffer; gravity helps buffer to pull chromatin down the slide. (Note: It usually takes me 15-20 seconds to pull a slide completely out of the coplin jar.)
7. Fix in 4% paraformaldehyde/formaldehyde diluted in PBS, PBS-Tween (0.05%) or PBST (0.1% Triton X-100) for 10 minutes at RT.
8. Rinse slides in 1X PBS+0.05% Tween-20 for 3 minutes at RT.

9. Block slides for 10-30 minutes at RT with PBS-Triton (1XPBS + 1%BSA + 0.5%Triton X-100 + 0.02% NaAzide) or your favorite blocking reagent.
10. Add primary antibody and incubate for 2-3 hours at RT or **overnight at 4°C**.
Overnight incubation at 4°C significantly decreases background; background needs to be minimized on the fiber preps since, ideally, the fibers should be very thin and antibody binding needs to be clearly distinguished from any background fluorescence.
11. Wash slides in PBS-Tween (1X PBS + 0.05%Tween), 3 times for 5 minutes. Do not shake coplin jar as that will dislodge fibers from the slide.
12. Add secondary antibody and incubate for 1 hour at RT. I recommend using Alexa-Fluor conjugated secondaries – they give very strong fluorescence with little background. Ex. Alexa-488 (in place of FITC), Alexa-555, 568, 594 (in place of Rhod-amine or Texas Red). Cy3 is also a great conjugate for secondaries.
13. Wash as in step 7.
14. Counterstain with antifade (Vectashield) containing DAPI. The DAPI concentration should be 2-5ug/mL in antifade so that you can see the fibers.
15. If proceeding to do FISH or thymidine analog detection, crosslink antibodies and proteins in 4-10% formalin in KCM or PBS, 15 minutes at RT.
16. Rinse slides in PBS-T or 2XSSC (if necessary, slides can be stored for a few days at 4C in these buffers before doing FISH).

Salt-Detergent-Urea Lysis Buffer - make fresh each (important!)

25 mM Tris, pH7.5

0.5M NaCl

1% Triton X-100

0.5M Urea (Dilute Urea from 1M Urea that is made up fresh each time in dH₂O)