

This file provides information about the dynamic image sets in the folders QIBA_DYN_2p_v10_beta1_GE and QIBA_DYN_2p_v10_beta1_Siemens.

This data set provides noise-free images similar to those of version 5, except that a magnetic field strength of 3T is simulated for both the dynamic and T1 mapping series.

The image data for the GE and Siemens sets are identical. The only difference is that the GE series contains DICOM headers that a GE machine would generate, while the Siemens series contains DICOM headers from a Siemens machine.

Two sets of images are provided for each series. DICOM part 10 format images are in the DICOM directory. XML files are in the XML directory. The XML images allow the values for the DICOM tags to be altered using a text editor, and new DICOM images can then be generated using dcm4che's tool called "xml2dcm," available at <http://www.dcm4che.org/confluence/display/d2/dcm4che2+DICOM+Toolkit>.

The parameters used to generate this data are:

Magnetic strength = 3 T
Flip angle = 25 degrees
Repetition Time = 5 msec
Time interval between the DCE images = 0.5 seconds
Assumed T1 (in tissue) = 1500 msec
Assumed equilibrium magnetization (in tissue) = 50000
Assumed T1 (in blood vessel) = 1932 msec
Assumed equilibrium magnetization (in blood vessel) = 50000

GE timing information is included in the DICOM headers at fields 0008,0032 (Acquisition Time) and 0018,1060 (Trigger Time). The imaging start time for Siemens machines is given in the DICOM headers at fields 0008,0030 (Study Time) and 0008,0031 (Series Time), and acquisition time for each image is given in the DICOM headers at fields 0008,0032 (Acquisition Time) and 0008,0033 (Image Time).

The input function was derived from the JSim model ToftsKermode_two_parameter_3T_fa25_20150112.proj (available on the QIBA page of our website). A link to the JSim website is also provided. This function representing a blood concentration time curve was converted to a plasma concentration time curve assuming a blood hematocrit of 45%. The relaxivity of the gadolinium contrast agent at 3 T was assumed to be $0.0037 \text{ mmol}^{-1} \text{ msec}^{-1}$ (Rohrer M, Bauer H, et. al. Comparison of Magnetic Properties of MRI Contrast Media Solutions at Different Magnetic Field Strengths. Investigative Radiology. 2005 Nov;40(11):715-24. PubMed PMID: 16230904).

The data in the test image is organized as follows:

The test data is generated using several combinations of Ktrans and Ve, using the modified Tofts Kermode 2-parameter model. The Ktrans takes values {0.01, 0.02, 0.05, 0.1, 0.2, 0.35}. The Ve takes values {0.01, 0.05, 0.1, 0.2, 0.5}.

The test data contains 10×10 pixels patches of each Ktrans and Ve combination. While generating the test data, the Ve values {0.01, 0.05, 0.1, 0.2, 0.5} vary along the x direction. Ktrans values {0.01, 0.02, 0.05, 0.1, 0.2, 0.35} vary along the y direction.

The Vascular region of interest is the bottom 50×10 pixels strip of the image. The peak of the Vascular region is the top-left 25×10 pixels strip of the image. This strip also contains time point labels given in seconds. The Zero patch (Ktrans=0.0, Ve=0.5) is the top-right 25×10 pixels strip of the image.

The following is a detailed list giving the specific Ktrans, Ve combination used to generate each 10*10 pixel patch. The x,y location specifies the upper-left corner of each 10*10 pixel patch containing a specific Ktrans, Ve combination.

x	y	Ktrans	Ve
0	10	0.01	0.01
0	20	0.02	0.01
0	30	0.05	0.01
0	40	0.1	0.01
0	50	0.2	0.01
0	60	0.35	0.01
10	10	0.01	0.05
10	20	0.02	0.05
10	30	0.05	0.05
10	40	0.1	0.05
10	50	0.2	0.05
10	60	0.35	0.05
20	10	0.01	0.1
20	20	0.02	0.1
20	30	0.05	0.1
20	40	0.1	0.1
20	50	0.2	0.1
20	60	0.35	0.1
30	10	0.01	0.2
30	20	0.02	0.2
30	30	0.05	0.2
30	40	0.1	0.2
30	50	0.2	0.2
30	60	0.35	0.2
40	10	0.01	0.5
40	20	0.02	0.5
40	30	0.05	0.5
40	40	0.1	0.5
40	50	0.2	0.5
40	60	0.35	0.5